


## Product Datasheet

### ATP5F1B antibody (orb1325096)

|                            |                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to ATP5F1B                                                       |
| <b>Species/Host</b>        | Mouse                                                                                      |
| <b>Reactivity</b>          | Canine, Human, Monkey, Rat                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                               |
| <b>Tested Applications</b> | FC, IF, IHC, WB                                                                            |
| <b>Immunogen</b>           | Full length human recombinant protein of human ATP5B (NP_001677) produced in HEK293T cell. |
| <b>Target</b>              | ATP5F1B                                                                                    |
| <b>Note</b>                | For research use only                                                                      |
| <b>Isotype</b>             | IgG1                                                                                       |
| <b>Clonality</b>           | Monoclonal                                                                                 |
| <b>MW</b>                  | 56.6 kDa                                                                                   |
| <b>Uniprot ID</b>          | <a href="#">P06576</a>                                                                     |
| <b>NCBI</b>                | <a href="#">NM_001686</a>                                                                  |
| <b>Dilution Range</b>      | WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                            |



Immunohistochemical staining of paraffin...