


Product Datasheet

MAN1 (LEMD3) antibody (orb1325135)

Description	Mouse monoclonal antibody to MAN1 (LEMD3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Full length human recombinant protein of human LEMD3 (NP_055134) produced in HEK293T cell.
Target	LEMD3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	100.0 kDa
Uniprot ID	Q9Y2U8
NCBI	NM_014319
Dilution Range	WB 1:1000, IF 1:100, Flow 1:100, IP 2-4ug/mg
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...