


Product Datasheet

ID3 antibody (orb1325140)

Description	Mouse monoclonal antibody to ID3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human ID3 (NP_002158) produced in HEK293T cell.
Target	ID3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	12.8 kDa
Uniprot ID	Q02535
NCBI	NM_002167
Dilution Range	WB 1:500, IHC 1:50, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...