


Product Datasheet

Sorbitol Dehydrogenase (SORD) antibody (orb1325146)

Description	Mouse monoclonal antibody to Sorbitol Dehydrogenase (SORD)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Full length human recombinant protein of human SORD (NP_003095) produced in HEK293T cell.
Target	SORD
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	38.3 kDa
Uniprot ID	Q00796
NCBI	NM_003104
Dilution Range	IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent staining of COS7 cells...