


Product Datasheet

Livin (BIRC7) antibody (orb1325169)

Description	Mouse monoclonal antibody to Livin (BIRC7)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human BIRC7 (NP_071444) produced in HEK293T cell.
Target	BIRC7
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	32.8 kDa
Uniprot ID	Q96CA5
NCBI	NM_139317
Dilution Range	WB 1:2000, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...