

Product Datasheet

Cytochrome P450 2E1 (CYP2E1) antibody (orb1325209)

Description Mouse monoclonal antibody to Cytochrome P450 2E1 (CYP2E1)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IF, IHC, WB

Immunogen Full length human recombinant protein of human CYP2E1 (NP_000764) produced in HEK293T cell.

Target CYP2E1

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 56.8 kDa

Uniprot ID [P05181](#)

NCBI [NM_000773](#)

Dilution Range WB 1:500, IHC 1:50, IF 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...