


Product Datasheet

IBRDC2 (RNF144B) antibody (orb1325227)

Description	Mouse monoclonal antibody to IBRDC2 (RNF144B)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-256 of human RNF144B (NP_877434) produced in E. coli.
Target	RNF144B
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	33.5 kDa
Uniprot ID	Q7Z419
NCBI	NM_182757
Dilution Range	WB 1:2000, IF 1:50~100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...