


Product Datasheet

SLC2A6 antibody (orb1325278)

Description	Mouse monoclonal antibody to SLC2A6
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human SLC2A6 expression vector
Target	SLC2A6
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	54.5 kDa
Uniprot ID	Q9UGQ3
NCBI	NM_017585
Dilution Range	IHC 1:200, WB: 1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...