


Product Datasheet

PP5 (PPP5C) antibody (orb1325296)

Description	Mouse monoclonal antibody to PP5 (PPP5C)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human PPP5C expression vector
Target	PPP5C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	56.9 kDa
Uniprot ID	P53041
NCBI	NM_006247
Dilution Range	WB 1:1000, IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...