


Product Datasheet

SSX2 antibody (orb1325297)

Description	Mouse monoclonal antibody to SSX2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human SSX2 expression vector
Target	SSX2
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.6 kDa
Uniprot ID	Q16385
NCBI	NM_003147
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...