


Product Datasheet

BubR1 (BUB1B) antibody (orb1325306)

Description	Mouse monoclonal antibody to BubR1 (BUB1B)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human BUB1B expression vector
Target	BUB1B
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	119.5 kDa
Uniprot ID	O60566
NCBI	NM_001211
Dilution Range	WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...