


Product Datasheet

PP5 (PPP5C) antibody (orb1325321)

Description	Mouse monoclonal antibody to PP5 (PPP5C)
Species/Host	Mouse
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Full length human recombinant protein of human PPP5C (NP_006238) produced in HEK293T cell.
Target	PPP5C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	56.7 kDa
Uniprot ID	P53041
NCBI	NM_006247
Dilution Range	WB 1:1000~2000, IHC 1:50, IF 1:50~100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...