

Product Datasheet

Cytochrome P450 Reductase (POR) antibody (orb1325342)

Description Mouse monoclonal antibody to Cytochrome P450 Reductase (POR)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications FC, IF, IHC, WB

Immunogen Full-length protein expressed in 293T cell transfected with human POR expression vector

Target POR

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 76.7 kDa

Uniprot ID [P16435](#)

NCBI [NM_000941](#)

Dilution Range WB 1:1000, IHC 1:50, IF 1:100, Flow 1:100

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...