

Product Datasheet

Protein kinase Y linked (PRKY) antibody (orb1325351)

Description Mouse monoclonal antibody to Protein kinase Y linked (PRKY)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IF

Immunogen Full length human recombinant protein of human PRKY (NP_002751) produced in HEK293T cell.

Target PRKY

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 31.7 kDa

Uniprot ID [O43930](#)

NCBI [NM_002760](#)

Dilution Range IF 1:100, Flow 1:100

Expiration Date 12 months from date of receipt.



immunofluorescent
staining of COS7
cells...