


Product Datasheet

DHFR antibody (orb1325371)

Description	Mouse monoclonal antibody to DHFR
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human DHFR expression vector
Target	DHFR
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.5 kDa
Uniprot ID	P00374
NCBI	NM_000791
Dilution Range	WB 1:2000 IHC 1:50 IF 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...