


Product Datasheet

LAT2 (SLC7A8) antibody (orb1325400)

Description	Mouse monoclonal antibody to LAT2 (SLC7A8)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human SLC7A8(NP_036376) produced in HEK293T cell.
Target	SLC7A8
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	58.2 kDa
Uniprot ID	Q9UHI5
NCBI	NM_012244
Dilution Range	WB 1:500~2000, IF 1:100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...