


Product Datasheet

Superoxide Dismutase 1 (SOD1) antibody (orb1325419)

Description	Mouse monoclonal antibody to Superoxide Dismutase 1 (SOD1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	Full length human recombinant protein of human SOD1 (NP_000445) produced in HEK293T cell.
Target	SOD1
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	15.8 kDa
Uniprot ID	P00441
NCBI	NM_000454
Dilution Range	WB 1:1000~2000, IHC 1:50, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...