


Product Datasheet

SCY1 like 3 (SCYL3) antibody (orb1325421)

Description	Mouse monoclonal antibody to SCY1 like 3 (SCYL3)
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human SCYL3 expression vector
Target	SCYL3
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	82.8 kDa
Uniprot ID	Q8IZE3
NCBI	NM_020423
Dilution Range	WB 1:2000 IHC 1:50 IF 1:100 Flow 1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...