


Product Datasheet

CD4 antibody (orb1325431)

Description	Mouse monoclonal antibody to CD4
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human CD4 (NP_000607) produced in HEK293T cell.
Target	CD4
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	51.1 kDa
Uniprot ID	P01730
NCBI	NM_000616
Dilution Range	WB 1:1000~2000, FC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...