


Product Datasheet

B Raf (BRAF) antibody (orb1325466)

Description	Mouse monoclonal antibody to B Raf (BRAF)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human BRAF expression vector
Target	BRAF
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	84.4 kDa
Uniprot ID	P15056
NCBI	NM_004333
Dilution Range	WB 1:2000 IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...