


Product Datasheet

SSB antibody (orb1325483)

Description	Mouse monoclonal antibody to SSB
Species/Host	Mouse
Reactivity	Canine, Human, Monkey
Conjugation	Unconjugated
Tested Applications	IF, IHC, IP, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human SSB expression vector
Target	SSB
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46.8 kDa
Uniprot ID	P05455
NCBI	NM_003142
Dilution Range	WB 1:1000~2000, IHC 1:150, IF 1:100, IP 2ug/500ul
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...