


Product Datasheet

STK3 antibody (orb1325499)

Description	Mouse monoclonal antibody to STK3
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human STK3 expression vector
Target	STK3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	56.3 kDa
Uniprot ID	Q13188
NCBI	NM_006281
Dilution Range	WB 1:2000; IHC 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...