


Product Datasheet

Chk2 (CHEK2) antibody (orb1325502)

Description	Mouse monoclonal antibody to Chk2 (CHEK2)
Species/Host	Mouse
Reactivity	Canine, Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Full length human recombinant protein of human CHEK2 (NP_009125) produced in HEK293T cell.
Target	CHEK2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	60.9 kDa
Uniprot ID	O96017
NCBI	NM_007194
Dilution Range	WB 1:500~1000, IF 1:50~100, FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
HEK293T
cells
transfected...