


Product Datasheet

GFAP antibody (orb1325558)

Description	Mouse monoclonal antibody to GFAP
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Full-length protein expressed in 293T cell transfected with human GFAP expression vector
Target	GFAP
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	49.9 kDa
Uniprot ID	P14136
NCBI	NM_002055
Dilution Range	WB 1:2000, IHC: 1:50-1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...