



Product Datasheet

DDK (FLAG) antibody (orb1325625)


Description	Mouse monoclonal antibody to DDK (FLAG)
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, IP, WB
Immunogen	Anti-DDK monoclonal antibody is produced by immunizing mice with a synthetic peptide coupled to KLH
Target	DDK (FLAG)
Preservatives	PBS (pH7.3) containing 1% rAlbumin, 50% glycerol and 0.02% sodium azide
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	WB: 1:2000, IP: 1:200, IF 1:1000, FC: 1:1000, IHC: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of immunoblot analysis...



Immunoblot analysis of randomly picked 7...



Western Blot analysis of DDK (FLAG) fuse...