


Product Datasheet

Neurogenin3 (NEUROG3) antibody (orb1325631)

Description	Mouse monoclonal antibody to Neurogenin3 (NEUROG3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant fragment expressed in E.coli corresponding to amino acids 1-214 of human Ngn3
Target	NEUROG3
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	22.9 kDa
Uniprot ID	Q9Y4Z2
NCBI	NM_020999
Dilution Range	WB: 1:500, IHC: 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...