


Product Datasheet

MAPK6 antibody (orb1325636)

Description	Mouse monoclonal antibody to MAPK6
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 345-721 of human ERK3
Target	MAPK6
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	82.5 kDa
Uniprot ID	Q16659
NCBI	NM_002748
Dilution Range	WB: 1:1000, IF (1:100)
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...