


Product Datasheet

EPCAM antibody (orb1325654)

Description	Mouse monoclonal antibody to EPCAM
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein expressed in E.coli corresponding to amino acids 24-265 of human EpCAM
Target	EPCAM
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	32.7 kDa
Uniprot ID	P16422
NCBI	NM_002354
Dilution Range	WB: 1:200 - 1:500, IF: 1:100
Expiration Date	12 months from date of receipt.



immunofluorescent
staining of HeLa
cells...