


Product Datasheet

Neurotrophin 4 (NTF4) antibody (orb1325675)

Description	Mouse monoclonal antibody to Neurotrophin 4 (NTF4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 81-210 human NT4 (NP_006170) produced in E.coli.
Target	NTF4
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	22.2 kDa
Uniprot ID	P34130
NCBI	NM_006179
Dilution Range	WB 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...