

Product Datasheet

GFAP antibody (orb1325733)

Description	Mouse monoclonal antibody to GFAP
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human recombinant protein fragment corresponding to amino acids 292-432 of human GFAP (NP_002046) produced in E.coli.
Target	GFAP
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	49.9 kDa
Uniprot ID	P14136
NCBI	NM_002055
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.