

Product Datasheet

Neurofilament (NEFL) antibody (orb1325742)

Description	Mouse monoclonal antibody to Neurofilament (NEFL)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human NEFL (NP_006149). The exact sequence is proprietary.
Target	NEFL
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	61.5 kDa
Uniprot ID	P07196
NCBI	NM_006158
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.