

Product Datasheet

IL17 (IL17A) antibody (orb1325763)

Description	Mouse monoclonal antibody to IL17 (IL17A)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full length human recombinant protein of human IL-17A (NP_002181) produced in 293 cell.
Target	IL17A
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	17.5 kDa
Uniprot ID	Q16552
NCBI	NM_002190
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.