


Product Datasheet

SARS-CoV-2 S Protein antibody (orb1325814)

Description	Mouse monoclonal antibody to SARS-CoV-2 S Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Spike-RBD protein produced in E.coli.
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	WB 1:500-1000, ELISA 1:5000-10001
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant SAR...