

Product Datasheet

SARS-CoV-2 N Protein antibody (orb1325827)

Description	Mouse monoclonal antibody to SARS-CoV-2 N Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full length recombinant SARS-Cov-2 N protein produced in E.coli.
Target	N Protein
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
MW	47kDa
NCBI	YP_009724397
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.