


Product Datasheet

SARS-CoV-2 N Protein antibody (orb1325839)

Description	Mouse monoclonal antibody to SARS-CoV-2 N Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant full length SARS-CoV-2 N protein produced in E.coli.
Target	N Protein
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	WB 1:500-1000, ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant full length SARS-CoV-2 N protein produced in E.coli.