

Product Datasheet

S100A8/A9 antibody (orb1325857)

Description	Mouse monoclonal antibody to S100A8/A9
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full length human S100A8/A9 co-expression produced in E.coli.
Target	S100A8/A9
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	25.76KD
Dilution Range	ELISA 1:1000-5000
Expiration Date	12 months from date of receipt.