

Product Datasheet

NGAL (LCN2) antibody (orb1325889)

Description	Mouse monoclonal antibody to NGAL (LCN2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human recombinant protein fragment corresponding to amino acids 21-198 of human NGAL (NP_005555) produced in E.coli.
Target	LCN2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	22.6 kDa
Uniprot ID	P80188
NCBI	NM_005564
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.