

Product Datasheet

IgG4 antibody (orb1325918)

Description	Mouse monoclonal antibody to IgG4
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Native Human IgG4 protein
Target	IgG4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	150KD
Dilution Range	ELISA 1:1000-5000
Expiration Date	12 months from date of receipt.