


Product Datasheet

CD21 (CR2) antibody (orb1325930)

Description	Mouse monoclonal antibody to CD21 (CR2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 208-587 of human CD21 (NP_001868) produced in E.coli.
Target	CR2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	113 kDa
Uniprot ID	P20023
NCBI	NM_001877
Dilution Range	IHC 1:250
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...