

Product Datasheet

D-Dimer antibody (orb1326083)

Description	Mouse monoclonal antibody to D-Dimer
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Native Human D-Dimer protein
Target	D-Dimer
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	180 kDa
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.