

Product Datasheet

Calponin (CNN1) antibody (orb1326104)

Description Mouse monoclonal antibody to Calponin (CNN1)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications IHC

Immunogen Human recombinant protein fragment corresponding to amino acids 67-257 of human Calponin1 (NP_001290) produced in E.coli.

Target CNN1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 33.2 kDa

Uniprot ID [P51911](#)

NCBI [NM_001299](#)

Dilution Range IHC 1:800

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...