


Product Datasheet

IL3RA antibody (orb1326144)

Description	Mouse monoclonal antibody to IL3RA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Full length human recombinant protein of human IL3RA (NP_002174) produced in HEK293T cell.
Target	IL3RA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	43.33 kDa
Uniprot ID	P26951
NCBI	NM_002183
Dilution Range	WB 1:250-1000, FLOW 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...