


Product Datasheet

Syndecan 1 (SDC1) antibody (orb1326149)

Description	Mouse monoclonal antibody to Syndecan 1 (SDC1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Human recombinant protein fragment corresponding to amino acids 23-254 of human SDC1 (NP_001006947) produced in E.coli.
Target	SDC1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	32.5 kDa
Uniprot ID	P18827
NCBI	NM_001006946
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of living 293T ...