

Product Datasheet

Syndecan 1 (SDC1) antibody (orb1326155)

Description Mouse monoclonal antibody to Syndecan 1 (SDC1)

Species/Host Mouse

Reactivity Human

Conjugation Unconjugated

Tested Applications IHC

Immunogen Full length human recombinant protein of human SDC1 (NP_001006947) produced in HEK293T cell.

Target SDC1

Note For research use only

Application notes For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.

Isotype IgG1

Clonality Monoclonal


MW 32.5 kDa

Uniprot ID [P18827](#)

NCBI [NM_001006946](#)

Dilution Range IHC 1:400

Expiration Date 12 months from date of receipt.



Immunohistochemical staining of paraffin...