


Product Datasheet

CD28 antibody (orb1326172)

Description	Mouse monoclonal antibody to CD28
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Full length human recombinant protein of human CD28 (NP_006130) produced in HEK293T cell.
Target	CD28
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	25.1 kDa
Uniprot ID	P10747
NCBI	NM_006139
Dilution Range	FLOW 1:100
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
living 293T ...