


Product Datasheet

iNOS (NOS2) antibody (orb1326201)

Description	Mouse monoclonal antibody to iNOS (NOS2)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 539-970 of human Nos2 (NP_000616) produced in E.coli.
Target	NOS2
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	130.9 kDa
Uniprot ID	P35228
NCBI	NM_000625
Dilution Range	WB 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...