


Product Datasheet

PKC alpha (PRKCA) antibody (orb1326218)

Description	Mouse monoclonal antibody to PKC alpha (PRKCA)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 350-672 of human PRKCA (NP_002728) produced in E.coli.
Target	PRKCA
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG2b
Clonality	Monoclonal
MW	76.6 kDa
Uniprot ID	P17252
NCBI	NM_002737
Dilution Range	WB 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...