


Product Datasheet

DDX4 antibody (orb1326228)

Description	Mouse monoclonal antibody to DDX4
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-400 of human DDX4 (NP_077726) produced in E.coli.
Target	DDX4
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	79.1 kDa
Uniprot ID	Q9NQI0
NCBI	NM_024415
Dilution Range	WB 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...