


Product Datasheet

IL10 antibody (orb1326252)

Description	Mouse monoclonal antibody to IL10
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 19-178 of human IL10 (NP_000563) produced in E.coli.
Target	IL10
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	20.52 kDa
Uniprot ID	P22301
NCBI	NM_000572
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



HEK293T cells were transfected with the ...