


Product Datasheet

ICOS antibody (orb1326318)

Description	Mouse monoclonal antibody to ICOS
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Full length human recombinant protein of human ICOS (NP_036224) produced in HEK293T cell.
Target	ICOS
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	20.2 kDa
Uniprot ID	Q9Y6W8
NCBI	NM_012092
Dilution Range	FLOW 1:100~400
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
living 293T ...